

NPY DATA FORMAT

NPY File Description for MLM-IV System

The `NPY` files generated by the [MLM-IV-SIMPLOT](#) program store simulated Langmuir probe I-V data in a format optimized for analysis. These files are specifically designed to be compatible with [PlasmaPy](#) and our proprietary [MLM-IV-ANALYSIS](#) program for detailed diagnostics.

File Structure

1. **File Format:** NPY [NumPy array format](#).

- The file contains a 2D NumPy array with two rows:
 - Row 1:** Probe voltage values (in volts, `V`).
 - Row 2:** Corresponding probe current values (in amperes, `A`).

2. **Purpose:**

- [MLM-IV-SIMPLOT](#): Simulates Langmuir probe I-V curves, including physical effects such as electron/ion leakage currents, noise, and plasma temperature. The results are saved in the `NPY` format for consistency and compatibility.
- [MLM-IV-ANALYSIS](#): Processes the generated I-V curves for plasma parameters such as electron temperature (`Te`), plasma potential (`Vp`), and electron density (`ne`), using methods like derivative analysis and curve fitting.

3. **Interoperability:**

- PlasmaPy can directly read and analyze these files, leveraging its plasma physics functions for extended studies.
- [MLM-IV-ANALYSIS](#) enhances the analysis by applying domain-specific diagnostics, such as identifying plasma potential from I-V derivatives or calculating electron densities.

Example Data Usage

- The first row provides the voltage range of the Langmuir probe experiment, typically spanning a defined range (e.g., -20V to 20V).
- The second row contains the corresponding measured or simulated current values, reflecting the plasma interaction with the probe under varying bias conditions.

Compatibility

- The `NPY` format's compact structure ensures that both PlasmaPy and custom analysis tools can efficiently load, process, and visualize the I-V data.
- This format supports rapid integration into workflows where raw data from [MLM-IV-SIMPLOT](#) is directly analyzed, fitted, and visualized for real-time insights or archival purposes.

Generated Files

- **Raw Simulated Data:** Contains the I-V curve with added noise and physical corrections.
- **Processed Data:** May include averaged noisy samples, theoretical curves, or adjusted fits, depending on the operation mode of the program.

Key Features

- **Efficient Storage:** Binary format ensures small file size and quick access.
- **Ready-to-Use Data:** Structured for direct ingestion into [PlasmaPy](#) or [MLM-IV-ANALYSIS](#).
- **Detailed Metadata:** Captures experimental parameters implicitly via the consistent voltage-current structure.

This format ensures robust compatibility for plasma diagnostics across various tools and platforms, enhancing the reproducibility and reliability of Langmuir probe experiments.